



Product Datasheet

JUND (phospho-S255) antibody (orb338848)

Description	Rabbit polyclonal antibody to JUND
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S255 of human JUND protein. The exact sequence is proprietary.
Target	JUND
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P15066, P17535, P52909
Entrez	24518, 16478, 3727
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of JUND (Phospho-S255) antibody.



Immunohistochemical analysis of JUND (Phospho-S255) antibody.