



## Product Datasheet

### PKC iota/zeta antibody (orb338857)


|                            |                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PRKCI                                                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rabbit, Rat, Zebrafish                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                        |
| <b>Tested Applications</b> | IF, IH, WB                                                                                                                                          |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human PKC iota/zeta. The exact sequence is proprietary.        |
| <b>Target</b>              | PRKCI; PRKCZ                                                                                                                                        |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                           |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                           |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                    |
| <b>Note</b>                | For research use only                                                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                                                                          |
| <b>Source</b>              | Rabbit                                                                                                                                              |
| <b>Uniprot ID</b>          | <a href="#">Q62074</a> , <a href="#">Q05513</a> , <a href="#">P09217</a> , <a href="#">Q02956</a> , <a href="#">F1M7Y5</a> , <a href="#">P41743</a> |
| <b>Entrez</b>              | <a href="#">5584</a> , <a href="#">84006</a> , <a href="#">18762</a> , <a href="#">18759</a> , <a href="#">5590</a> , <a href="#">25522</a>         |
| <b>Dilution Range</b>      | WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                     |



Western blot analysis of PKC iota/zeta e...



Immunohistochemical analysis of PKC iota...



Immunofluorescent analysis of PKC iota/z...