



Product Datasheet

p38 (phospho-Y182) antibody (orb338863)


Description	Rabbit polyclonal antibody to MAPK14
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y182 of human p38 protein. The exact sequence is proprietary.
Target	MAPK14
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q16539 , P47811 , P70618
Entrez	1432 , 26416
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of p38 (Phospho-Y182) antibody.



Immunohistochemical analysis of p38 (Phospho-Y182) antibody.



Immunofluorescent analysis of p38 (Phospho-Y182) antibody.