



Product Datasheet

RB1 (phospho-S807) antibody (orb338869)


Description	Rabbit polyclonal antibody to RB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S807 of human RB1 protein. The exact sequence is proprietary.
Target	RB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P13405, P33568, P06400
Entrez	24708, 5925, 19645
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of RB1 (Phospho-S807) antibody.



Immunohistochemical analysis of RB1 (Phospho-S807) antibody.



Immunofluorescent analysis of RB1 (Phospho-S807) antibody.