



Product Datasheet

EPHB1/2 antibody (orb338937)

Description	Rabbit polyclonal antibody to EPHB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human EPHB1/2. The exact sequence is proprietary.
Target	EPHB1; EPHB2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P29323 , P54762 , Q8CBF3 , P09759 , P54763
Entrez	24338 , 270190 , 2047 , 1969
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.




Western blot analysis of EPHB1/2 express...



Immunohistochemical analysis of EPHB1/2

...



Immunofluorescent analysis of EPHB1/2 st...