

Product Datasheet

Histone H3 (phospho-T3) antibody (orb338957)

Description

Rabbit polyclonal antibody to HIST1H3A

Species/Host

Rabbit

Reactivity


Human, Mouse, Porcine, Primate, Rabbit, Rat, Zebrafish

Conjugation

Unconjugated

Tested

IH, WB

Applications


Western blot analysis of Histone H3 (Pho...)

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T3 of human Histone H3 protein. The exact sequence is proprietary.

Target

HIST1H3A

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source


Rabbit

Uniprot ID
P68431, P68433, Q6LED0
Entrez
8352, 8356, 8354, 8968, 319153, 8357, 97908, 682330, 319152, 8358, 8355, 360198, 291159, 8350, 8353, 8351
Dilution Range

WB: 1:500-1000, IHC-P: 1:100-200

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of Histone H3