




Product Datasheet

IKB alpha antibody (orb33904)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NFKBIA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35609
Uniprot ID	P25963
NCBI	NP_065390.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis of THP-1 whole cell lysate L...



Western blot analysis of Daudi whole cel...