


Product Datasheet

XPO1 (AcK568) antibody (orb339060)

Description	Rabbit polyclonal antibody to XPO1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K568 of human XPO1 protein. The exact sequence is proprietary.
Target	XPO1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q80U96, Q6P5F9, O14980
Entrez	103573, 85252, 7514
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of XPO1 (Acetyl-K568) antibody.