


Product Datasheet

14-3-3-pan (AcK51/49) antibody (orb339061)

Description	Rabbit polyclonal antibody to YWHAB
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K51/49 of human 14-3-3-pan protein. The exact sequence is proprietary.
Target	YWHAB; YWHAG; YWHAQ; YWHAZ; SFN
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P61983, P27348, P31947, P63104, P35213, P63101, Q9CQV8, P68254, P31946, P68255, P61981, P63102, P61982
Entrez	7529, 22630, 7534, 22631, 10971, 22628, 7532, 25577, 2810, 54401, 56011, 56010, 25578
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 14-3-3-pan (AcK51/49) antibody.