

Product Datasheet

Alpha-tubulin (AcK352) antibody (orb339083)

Description

Rabbit polyclonal antibody to TUBA1A

Species/Host

Rabbit

Reactivity

Human, Mouse, Primate, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K352 of human Alpha-tubulin protein. The exact sequence is proprietary.

Target

TUBA1A; TUBA1B; TUBA1C; TUBA3C; TUBA4A; TUBA8

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P68369](#), [Q71U36](#), [P68366](#), [P68370](#), [Q13748](#), [Q9BQE3](#), [Q6AYZ1](#), [Q6P9V9](#), [P05213](#), [P68373](#), [Q9NY65](#), [P68363](#)

Entrez

[113457](#), [51807](#), [22146](#), [22142](#), [500929](#), [64158](#), [7278](#), [100909441](#), [1962](#), [22143](#), [7277](#), [84790](#), [10376](#), [7846](#)

Dilution Range


WB: 1:500-1000

Expiration Date

12 months from date of receipt.

kDa

A B



Western blot analysis of Alpha-tubulin (...)