

Product Datasheet

p53 antibody (orb339086)

Description	Rabbit polyclonal antibody to TP53
Species/Host	Rabbit
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic Mono-methylated peptide corresponding to residues surrounding K370 of human p53 protein. The exact sequence is proprietary.
Target	TP53
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P04637
Entrez	7157
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.

