



Product Datasheet

MC5 Receptor antibody (orb339174)

Description	Rabbit polyclonal antibody to MC5R
Species/Host	Rabbit
Reactivity	Human, Mouse, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human MC5 Receptor. The exact sequence is proprietary.
Target	MC5R
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P41149 , P33032
Entrez	4161
Dilution Range	WB: 1:500-1000, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of MC5 Receptor ex...



Immunofluorescent analysis of MC5 Recept...