



Product Datasheet

ATG4D antibody (orb33920)


Description	ATG4D antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Cysteine protease ATG4D protein (14-43AA)
Target	ATG4D
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86TL0
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(lane 1), H...



Immunohistochemical staining of human sm...



Immunofluorescence analysis of PC-3 cell...