

Product Datasheet

Factor XII HC peptide (orb339214)

Description	The peptide is used to Factor XII HC antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	F12
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338820 (Factor XII HC antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q80YC5, P00748
Entrez	2161, 58992
Expiration Date	12 months from date of receipt.