

Product Datasheet

HSD3B2 peptide (orb339215)

Description	The peptide is used to HSD3B2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HSD3B2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338821 (HSD3B2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26149, P22072, P26439
Entrez	682974, 3284, 15493
Expiration Date	12 months from date of receipt.