

Product Datasheet

GLUR4 peptide (orb339230)

Description	The peptide is used to GLUR4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	GRIA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338836 (GLUR4 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19493, P48058, Q9Z2W8
Entrez	29629, 2893, 14802
Expiration Date	12 months from date of receipt.