

Product Datasheet

mPR alpha peptide (orb339237)

Description	The peptide is used to mPR alpha antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PAQR7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338843 (mPR alpha antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q86WK9, Q80ZE4
Entrez	71904, 164091
Expiration Date	12 months from date of receipt.