

Product Datasheet

mPR gamma peptide (orb339238)

Description	The peptide is used to mPR gamma antibody reactivity.
Reactivity	Human, Mouse, Sheep
Tested Applications	Blocking
Target	PAQR5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338844 (mPR gamma antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NXK6, Q9DCU0
Entrez	74090, 54852
Expiration Date	12 months from date of receipt.