

Product Datasheet

JUND (phospho-S255) peptide (orb339242)

Description	The peptide is used to JUND (pS255) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	JUND
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338848 (JUND (phospho-S255) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15066, P17535, P52909
Entrez	16478, 3727, 24518
Expiration Date	12 months from date of receipt.