

Product Datasheet

p300 peptide (orb339244)

Description	The peptide is used to p300 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	EP300
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338850 (p300 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	B2RWS6, Q09472
Entrez	2033, 328572
Expiration Date	12 months from date of receipt.