

Product Datasheet

Ezrin/Radixin/Moesin peptide (orb339246)

Description	The peptide is used to Ezrin/Radixin/Moesin antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	MSN; RDX; EZR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338852 (Ezrin/Radixin/Moesin antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26038, P35241, P15311, P31977, P26043, O35763, P26040, P26041
Entrez	19684, 4478, 17698, 81521, 5962, 54319, 22350
Expiration Date	12 months from date of receipt.