

Product Datasheet

HNF4 alpha/gamma peptide (orb339255)

Description	The peptide is used to HNF4 alpha/gamma antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HNF4A; HNF4G
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338861 (HNF4 alpha/gamma antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WUU6 , P41235 , P22449 , Q14541 , P49698
Entrez	3172 , 30942 , 15378 , 25735 , 3174
Expiration Date	12 months from date of receipt.