

Product Datasheet

Histone H3 peptide (orb339256)

Description	The peptide is used to Histone H3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	HIST1H3A; HIST1H3B; HIST1H3C; HIST1H3D; HIST1H3E; HIST1H3F; HIST1H3G; HIST1H3H; HIST1H3I; HIST1H3J; HIST2H3A; HIST2H3C; HIST2H3D; H3F3A; H3F3B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338862 (Histone H3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68431, P84244, P84243, P84245, P68433, Q71DI3, Q6LED0, P84228
Entrez	15078, 333932, 8352, 8351, 291159, 8350, 8353, 3021, 126961, 8354, 8356, 653604, 15077, 8358, 8357, 8355, 8968, 100361558, 3020, 319152
Expiration Date	12 months from date of receipt.