

Product Datasheet

Estrogen Receptor alpha peptide (orb339260)

Description	The peptide is used to Estrogen Receptor alpha antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	ESR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338866 (Estrogen Receptor alpha antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P03372
Entrez	2099
Expiration Date	12 months from date of receipt.