

Product Datasheet

Syntaxin 1A (phospho-S14) peptide (orb339268)

Description	The peptide is used to Syntaxin 1A (pS14) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	STX1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338874 (Syntaxin 1A (phospho-S14) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16623, O35526, P32851
Entrez	116470, 20907, 6804
Expiration Date	12 months from date of receipt.