

Product Datasheet

SMAD2 peptide (orb339269)

Description	The peptide is used to SMAD2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	SMAD2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338875 (SMAD2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O70436, Q15796, Q62432
Entrez	17126, 4087, 29357
Expiration Date	12 months from date of receipt.