

Product Datasheet

MYPT1 peptide (orb339277)

Description	The peptide is used to MYPT1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PPP1R12A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338883 (MYPT1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9DBR7, Q10728, O14974
Entrez	17931, 116670, 4659
Expiration Date	12 months from date of receipt.