

Product Datasheet

PERK peptide (orb339280)

Description	The peptide is used to PERK antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	EIF2AK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338886 (PERK antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z2B5 , Q9NZJ5 , Q9Z1Z1
Entrez	29702 , 9451
Expiration Date	12 months from date of receipt.