

Product Datasheet

DUSP4 peptide (orb339285)

Description	The peptide is used to DUSP4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	DUSP4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338891 (DUSP4 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13115, Q62767, Q8BFV3
Entrez	60587, 319520, 1846
Expiration Date	12 months from date of receipt.