

Product Datasheet

DCP1A peptide (orb339287)

Description	The peptide is used to DCP1A antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	DCP1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338893 (DCP1A antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NPI6, Q91YD3
Entrez	75901, 55802
Expiration Date	12 months from date of receipt.