

Product Datasheet

MENA peptide (orb339291)

Description	The peptide is used to MENA antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ENAH
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338897 (MENA antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q03173, Q8N8S7
Entrez	13800, 55740
Expiration Date	12 months from date of receipt.