

Product Datasheet

Histone H4 (AcK5) peptide (orb339293)

Description	The peptide is used to Histone H4 (AcK5) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338899 (Histone H4 (AcK5) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62806, P62804, P62805
Entrez	319161, 320332, 100041230, 326620, 8368, 102557184, 8361, 680097, 554313, 8294, 8365, 500351, 102548682, 97122, 8364, 100912290, 502913, 8359, 8370, 100912418, 102641229, 291152, 319160, 8366, 100360950, 295277, 326619, 64627, 121504, 319158, 319159, 8363, 69386, 102551184, 8362, 100912564, 319155, 8367, 319157, 8360, 319156
Expiration Date	12 months from date of receipt.