

Product Datasheet

STMN1 (phospho-S38) peptide (orb339297)

Description	The peptide is used to STMN1 (pS38) antibody reactivity.
Reactivity	Human, Porcine, Primate, Rat
Tested Applications	Blocking
Target	STMN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338903 (STMN1 (phospho-S38) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P13668, P16949
Entrez	3925, 29332
Expiration Date	12 months from date of receipt.