

Product Datasheet

MET peptide (orb339298)

Description	The peptide is used to MET antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	MET
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338904 (c-Met antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08581, P97523, P16056
Entrez	24553, 4233
Expiration Date	12 months from date of receipt.