

## Product Datasheet

### PKC delta (phospho-S645) peptide (orb339302)

**Description** The peptide is used to PKC delta (pS645) antibody reactivity.

**Reactivity** Human, Primate

**Tested Applications** Blocking

**Target** PRKCD

**Form/Appearance** Lyophilized powder

**Storage** Shipped at 4°C. Store at -20°C for one year.

**Note** For research use only

**Application notes** This peptide was used as immunogen to produce orb338908 (PKC delta (phospho-S645) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.

**Purity** > 85%

**Source** Synthetic

**Uniprot ID** **Q05655**

**Entrez** **5580**

**Expiration Date** 12 months from date of receipt.