

Product Datasheet

FANCD2 peptide (orb339304)

Description	The peptide is used to FANCD2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FANCD2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338910 (FANCD2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q80V62, Q6IV68, Q9BXW9
Entrez	312641, 2177, 211651
Expiration Date	12 months from date of receipt.