

Product Datasheet

Alpha-1A Adrenergic Receptor peptide (orb339313)

Description	The peptide is used to Alpha-1A Adrenergic Receptor antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADRA1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338919 (Alpha-1A Adrenergic Receptor antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97718, P43140, P35348
Entrez	148, 29412, 11549
Expiration Date	12 months from date of receipt.