

Product Datasheet

c-Jun peptide (orb339320)

Description	The peptide is used to c-Jun antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	JUN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338926 (c-Jun antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05627, P17325, P05412
Entrez	24516, 3725, 16476
Expiration Date	12 months from date of receipt.