

Product Datasheet

AMPK alpha 1 peptide (orb339326)

Description	The peptide is used to AMPK alpha 1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PRKAA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338932 (AMPK alpha 1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5EG47 , P54645 , Q13131
Entrez	105787 , 5562 , 65248
Expiration Date	12 months from date of receipt.