

Product Datasheet

Beta-adducin peptide (orb339336)

Description	The peptide is used to Beta-adducin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADD2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338942 (Beta-adducin antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q05764, Q9QYB8, P35612
Entrez	24171, 119, 11519
Expiration Date	12 months from date of receipt.