

Product Datasheet

Neuroglycan C peptide (orb339338)

Description	The peptide is used to Neuroglycan C antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CSPG5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338944 (Neuroglycan C antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O95196, Q71M36, Q9ERQ6
Entrez	10675, 29873, 50568
Expiration Date	12 months from date of receipt.