

Product Datasheet

Neuropilin 2 peptide (orb339345)

Description	The peptide is used to Neuropilin 2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NRP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338951 (Neuropilin 2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O35276, O35375, O60462
Entrez	8828, 81527, 18187
Expiration Date	12 months from date of receipt.