

Product Datasheet

SHCA (phospho-Y427) peptide (orb339349)

Description	The peptide is used to SHCA (pY427) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SHC1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338955 (SHCA (phospho-Y427) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29353, Q5M824, P98083
Entrez	20416, 85385, 6464
Expiration Date	12 months from date of receipt.