

Product Datasheet

CCR5 (phospho-S349) peptide (orb339350)

Description	The peptide is used to CCR5 (pS349) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit
Tested Applications	Blocking
Target	CCR5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338956 (CCR5 (phospho-S349) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51682, P51681
Entrez	12774, 1234
Expiration Date	12 months from date of receipt.