

Product Datasheet

Histone H3 (phospho-T3) peptide (orb339351)

Description	The peptide is used to Histone H3 (pT3) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	HIST1H3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338957 (Histone H3 (phospho-T3) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6LED0, P68431, P68433
Entrez	97908, 8358, 8352, 8351, 8355, 8357, 8968, 291159, 682330, 8350, 8353, 8354, 8356, 319152, 360198, 319153
Expiration Date	12 months from date of receipt.