

Product Datasheet

PI3K p85 alpha (phospho-Y607) peptide (orb339359)

Description	The peptide is used to PI3K p85 alpha (pY607) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PIK3R1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338965 (PI3K p85 alpha (phospho-Y607) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P27986, P26450, Q63787
Entrez	18708, 5295, 25513
Expiration Date	12 months from date of receipt.