

Product Datasheet

BCL2 peptide (orb339362)

Description	BCL2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BCL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338968 (BCL2 (phospho-S70) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49950, P10417
Entrez	12043, 24224
Expiration Date	12 months from date of receipt.