

Product Datasheet

Integrin beta 1 peptide (orb339365)

Description	Integrin beta 1 peptide
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	ITGB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338971 (CD29 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49134, P09055, P05556
Entrez	24511, 3688, 16412
Expiration Date	12 months from date of receipt.