

Product Datasheet

PPP2CA peptide (orb339370)

Description	The peptide is used to PPP2CA antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	PPP2CA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338976 (PPP2CA antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P63330, P63331, P67775
Entrez	5515, 24672, 19052
Expiration Date	12 months from date of receipt.