

Product Datasheet

CD32b peptide (orb339372)

Description	CD32b peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	FCGR2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338978 (CD32b antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31994
Entrez	2213
Expiration Date	12 months from date of receipt.